

Outline

Aim 1


- 1.1 Characterize and correct for B1 field inhomogeneity
- 1.2 Reduce distortion due to B0 field inhomogeneity
- 1.3 Real-time correction of head motion during scanning
- 1.4 Characterize and correct for gradient distortions
- 1.5 Continue with the calibration of T1 acquisitions

Outline

Aim 1

- 1.1 Characterize and correct for B1 field inhomogeneity
- 1.2 Reduce distortion due to B0 field inhomogeneity
- 1.3 Real-time correction of head motion during scanning
- 1.4 Characterize and correct for gradient distortions
- 1.5 Continue with the calibration of T1 acquisitions


B_1 Inhomogeneity Dielectric Resonance at 3T



MP_RAGE




T_1 -weighted GRE
TR=20ms, FA=30deg




PD-weighted GRE
TR=20ms, FA=5deg

B_1 Inhomogeneity Dielectric Resonance Effects at 3T




Estimated B_1 Transmit Map at 3T

B₁ Inhomogeneity Correction for Dielectric Resonance Effects at 3T




MP-RAGE




B₁ corrected
Synthetic T1-weighted

B_1 Inhomogeneity

Correction for Dielectric Resonance Effects at 3T



B_1 corrected
Synthetic T_1 -weighted



B_1 corrected
Synthetic "PD-weighted"

Outline

Aim 1

- 1.1 Characterize and correct for B1 field inhomogeneity
- 1.2 Reduce distortion due to B0 field inhomogeneity
- 1.3 Real-time correction of head motion during scanning
- 1.4 Characterize and correct for gradient distortions
- 1.5 Continue with the calibration of T1 acquisitions

Spatial Distortion due to Magnetic Susceptibility (B_0) Pulse-Sequence Dependence (S/I readout direction)




Spatial Distortion due to Magnetic Susceptibility (B_0) Pulse-Sequence Dependence (I/S readout direction)






Multi-Echo 3-D FLASH/SPGR Sequence

5 deg
Flip angle




30 deg
Flip angle



Reduction of Scan-Rescan Variance


EYE MOTION

FLASH 30




FLOW

FLASH 30



Multi-echo FLASH 30




Multi-echo FLASH 30



Overlay threshold: 8%

Image intensity histogram
of relative test-retest variability




Outline


Aim 1

- 1.1 Characterize and correct for B1 field inhomogeneity
- 1.2 Reduce distortion due to B0 field inhomogeneity
- 1.3 Real-time correction of head motion during scanning
- 1.4 Characterize and correct for gradient distortions
- 1.5 Continue with the calibration of T1 acquisitions

Rigid body motion is calculated every TR by comparing the new clover-leaf navigator with a pre-measured map of k-space in the vicinity of the target navigator




Clover-leaf navigator path in k-space




K-space values in vicinity of octant navigator

Real-Time Motion Tracking using Cloverleaf Navigators





Without motion correction



With real-time motion correction

With navigator motion correction it was possible to collect high resolution structural scans of the amnesic patient H.M.



3D FLASH (TR=20 ms, TE=10 ms, 1.3 x 1 x 1.3 mm, Tacq=7:45, BW=160 Hz/pixel, 6 averages)

Outline

Aim 1


- 1.1 Characterize and correct for B1 field inhomogeneity
- 1.2 Reduce distortion due to B0 field inhomogeneity
- 1.3 Real-time correction of head motion during scanning
- 1.4 Characterize and correct for gradient distortions
- 1.5 Continue with the calibration of T1 acquisitions

Gradient Nonlinearity Distortion

Uncorrected Phantom Data




**GE Whole-Body
CRM NVi/CVi**




**Siemens Whole-Body
Symphony/Sonata**

Gradient Nonlinearity Distortion Corrected Phantom Data




**GE Whole-Body
CRM NVi/CVi**



**Siemens Whole-Body
Symphony/Sonata**

Gradient Nonlinearity Distortion Uncorrected



Gradient Nonlinearity Correction Corrected

