



## Mission

The Surgical Navigation and Robotics (SNR) Laboratory enables more effective and less invasive treatment in image-guided therapy.

Brigham and Women's Hospital, Boston      Nobuhiko Hata

## SNR Mission

- We fulfill our mission through a commitment to:
  - Inventing disruptive computer and engineering methods
  - Applying the developed technologies in actual clinical cases and delivering unique feedback to the scientific research community
  - Sharing the research data, software, and device design with industry and academic peers
  - Applying synergistic coupling to scientific disciplines unaware of or presently disconnected from image-guided therapy

Brigham and Women's Hospital, Boston      Nobuhiko Hata

## SNR Value

**Quick cycle**

```

graph TD
    Science[Science] <--> Application[Application]
    Application <--> Engineering[Engineering]
    Engineering <--> Science
          
```

**Engineers in OR**

Brigham and Women's Hospital, Boston      Nobuhiko Hata

## Team Science

```

graph TD
    NCIGT[NCIGT] <--> SNR[SNR Lab]
    Internal[Internal collaborators] <--> SNR
    External[External collaborators] <--> SNR
    NA-MIC[NA-MIC] <--> SNR
          
```

Brigham and Women's Hospital, Boston      Nobuhiko Hata

## Acknowledgements

- Drs. Jolesz, Tempany, Kikinis, Wells, Golby (Brigham and Women's Hospital, Boston)
- Dr. Fichtinger (Queens Univ.), Dr. Morikawa (Shiga Univ., Japan)
- Slicer ([www.slicer.org](http://www.slicer.org))
- NIH U41 U41RR019703 (PI: Jolesz, BWH)  
National Center for Image Guided Therapy
- NIH 5U54EB005149  
National Alliance for Medical Image Computing (PI: Kikinis, BWH)
- NIH 5P01CA067165 (PI: Jolesz, BWH)  
MRI-guided therapy
- NIH BRP R01-CA111288 (PI: Tempany, BWH)  
Enabling Technologies for MRI-Guided Prostate Interventions
- NSF EEC 9731748 (PI: Taylor, JHU) ERC Center for Integrated Surgical Systems Technology
- CIMIT, Boston, MA
- Swimming Capsule Endoscope
- Intelligent Surgical Instruments Project of METI, Japan

Brigham and Women's Hospital, Boston      Nobuhiko Hata

## Outline of talk

- Part 1: MRI Guided Therapy
- Part 2: Robot as precision targeting device
- Part 3: Future: Robot as enabling tool for new therapy options (moving organ, molecular image-guided therapy, etc.)
- Part 4: Engineering and resource sharing

Brigham and Women's Hospital, Boston

Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MRI-guided therapy


State of art  
Unmet needs

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Image Guided Therapy

- Guides procedures with re-sectioning of **pre-operative images** (MRI, CT, ...)
- Aims at
  - Fast and accurate
  - Surgery with minimized intervention
  - Better clinical outcome




Slide courtesy of Dr. Jolesz

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Navigation

- To see beyond skin
- Computer-assisted stereotactic neurosurgery [Kelly86]
- **Frameless stereotaxy** [Watanabe87]
  - With free-hand tracking device
- An **interactive** tool to correlate “pre-operative images” and “therapeutic action”




Subject for demonstration purpose only

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Issues in pre-operative image-based surgical navigation

- **Deformation** (breaks accurate correlation between image to intra-operative anatomy)
- Cannot capture the “effect” of the therapy



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Intra-operative MRI




Imaging modality to navigate, monitor, control surgery  
Slide courtesy of Dr. Jolesz Brigham and Women's Hospital, Boston, MA

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Intra-operative MRI

- MRI for
  - Planning
  - Targeting
  - Navigation
  - Monitoring




Dr. Morikawa

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Intra-operative MRI




pictures from the vendors' web site

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Navigation in MRI-guided therapy

- "Updated navigation"
  - Requires patient-to-image registration after each patient transport
- [Nimsky04] Philips open-bore 0.2 and close-bore 3T and Stealth Station
- BrainSUITE iMRI (Brain Lab/Siemens)




Philips

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MR-guided Liver ablation Therapy




Morikawa S, Hata et. Al., Acad Radiol. 2003;10(2):180-8  
Morikawa S, Hata et. Al., Acad Radiol. 2003;10(12):1442-9.

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MRI-guided Neurosurgery




Dimairo SP, Hata N, et al. IEEE Eng Med Biol Mag. 2006;25(5):67-73.

Brigham and Women's Hospital, Boston


Nobuhiko Hata

## fMRI and intraoperative MRI




Brigham and Women's Hospital, Boston

Nobuhiko Hata



## MR-guided Prostate Biopsy



Intra-procedural 0.5T FGR  
for needle control

Registered Pre-procedural  
0.5T T2w


Hata N et al, Radiology. 2001;220(1):263-8.

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Needle Placement Accuracy in Prostate Therapies

- Improved cancer detection in targeted MRI-guided biopsy [So 2005]
- Retrospective clinical analysis (N=10) showed 6.5mm error [Blumenfeld 2007]
- Need for better needle placement control
  - Understanding of needle-tissue interaction for trajectory planning, optimization and control



Blumenfeld P, Hata N, et al. J Magn Reson Imaging. 2007  
Sep;26(3):688-94.

Brigham and Women's Hospital, Boston

Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata


## MRI-compatible robot

Autonomous approach  
Precision Targeting

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Precision needle placement



- Prostate needle placement up to 6 mm inaccuracy
  - fixed needle guide template with holes spaced at least 5mm apart
  - limits needle trajectory position and orientation
- Robot as dynamic needle guide for accurate needle placement [Chinzei MICCAI 1999]

•AIST-Japan and Brigham and Women's Hospital

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MRI Robot



© MEL, AIST/jpn & SPL, BWH/I

- 5DOF translational stage
- Maneuvering two arms
- Accurate needle placement

DiMaio SP, Hata N, et. al. Comput Aided Surg. 2007;12(1):15-24.

Brigham and Women's Hospital, Boston

Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MRI compatible robot


- Comprehensive review by Elhawary et al. 2006
- Open bore robot
  - Prototype for prostate biopsy [Chinzei 1999] [DiMaio 2006]
  - Heart [Tajima 2003]
  - Liver ablation therapy [Hata 2008]
- Close bore robot
  - Brain [Masamune 1995]
  - Breast [Larson2004]
  - Prostate biopsy and brachytherapy [Fischer 2008]

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MRI compatible robot

- Precise needle placement by



•Dohi T, Hata N, et. al. Robotics in Computer Aided Surgery. Journal of Computer Aided Surgery. 1995;1(1):4-10

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## MR Compatibility

- Actuators: **Ultrasonic**, Pneumatic, Hydraulic motors
- Ceramics bearing
- Aluminum and Stainless
- Optical encoder: MRI compatible or custom made with detached optical pickup and circuit

Brigham and Women's Hospital, Boston


Nobuhiko Hata

## Our Proposal Needle-guiding robot

Steps

1. Volume imaging
2. Target definition
3. Remote Center of Motion (RCM) set on the target
4. Optimal path search

- Synergistic robots [Troccaz, Robotics Workshop CARS 03]
- Cooperative robots [Kumar, Taylor ICRA 2000, MICCAI03]




• Proof of concept by the robot and a retrospective liver case

*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*

## Synergistic Robot in MRI


- hands-on
  - Physician takes initiative
  - Robot assists targeting by keeping motion tip at the target
- Conceptually, integration of
  - Stereotactic device
  - Master-is-slave




*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*


## Configuration

Active 3DOF XYZ stage




Passive 2DOF Needle holder





*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*

## MRI Safe/Compatible




- Actuators
  - Ultrasonic motors
- Motion parts
  - Rigidity!
  - Stainless (SUS304) for screw rod and linear guide
  - Ceramics bearing
- Frame
  - Aluminum, Stainless
- Optical rotary and linear encoders
  - With counter circuit
  - Extension fiber and separated circuit

*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*


## Clinical applications!

- Concept Needle guiding robot [Hata et al. ICRA 2005] [Hata et al. JMIRI 2008]
- Three generation (in-lab machining -> professional robot re-design/assembly)
- August 2007: First clinical case in liver ablation therapy
- N=4 as of today





*J Magn Reson Imaging, 2008 Apr 11;27(5):1130-1138.*

*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*



*Brigham and Women's Hospital, Boston* *Nobuhiko Hata*




### Long-term Clinical Goal at BWH

- Diagnosis and therapy in multi-modality image setting, incl. high-field closed-bore scanner (3T)
  - High-quality imaging,
  - More prevalent in clinics and hospitals.
- Challenges
  - Access to lesions
  - Motion compensation
  - Device-tissue interaction

Brigham and Women's Hospital, Boston

Nobuhiko Hata



### In-bore Needle Placement

GE 3T Excite bore

Optical Encoder

Pneumatic Cylinder

Vertical Motion Mechanisms

Joint Project with Dr. Gabor Fichtinger, Johns Hopkins Univ.


Brigham and Women's Hospital, Boston

Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Vision Image Driven Robotics

Brigham and Women's Hospital, Boston

Nobuhiko Hata




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Moving organ: rationale


- Combining MRI and Robot
- Value
  - MRI can track motion of moving organ
  - Robot can synchronize its motion to the moving organ
- IGT of moving organ
  - Liver, Lung, Heart



•Ishikawa, Univ. Tokyo

Brigham and Women's Hospital, Boston

Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata



## Synergistic robot with motion compensation in MRI




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Navigator echo and projection profile matching

Tokuda, Hata, et.al. Acad Rad 2005




Brigham and Women's Hospital, Boston

x

Nobuhiko Hata

## Pulse Sequence

- FGRE + Navigator
  - Projections along two axes
  - 128 navigators / image
  - Matrix: 256x128 TR: 27ms~
  - Tracking time resolution = 4 x TR




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## System Integration to Signa/SP


- Off-the-shelf system + 3D Slicer



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Tumor Treatment of moving organ with FUS



InSightec

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Actively swimming capsule endoscope using MRI for energy delivery, imaging, and navigation

- Goal of CIMIT project
  - Novel swimming mechanism to effectively steer and position the capsule inside the small intestine
  - using MRI's magnetic field to wirelessly send energy and generate propulsion force
- Our multi-year clinical goal
  - To mount a miniature endoscope camera and treatment device onto our swimming capsule to aid in diagnosis and treatment of the small intestine.

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Benefit


- The capsule endoscope today is purely diagnostic
  - Cannot be used to take biopsies, apply therapy, or mark abnormalities for surgery.
  - Cannot be controlled once it has been ingested
  - Once it has passed a suspicious abnormality, its progress cannot be slowed to better visualize the area



- The proposed system allows the gastroenterologist to maneuver it using MRI along with guidance and monitoring mapping

Brigham and Women's Hospital, Boston

Nobuhiko Hata




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Background

Flagellar movement of microorganisms in low Reynolds number hydrodynamic field [Colgate et al]  
 the oscillating beam can create approximated sinusoidal traveling wave in viscous flow and produces propulsion force effectively (Kosa et al)



\*Colgate, J.E. and Lynch, K.M., IEEE Journal of Oceanic Engineering 29, (2004) 660-673.  
 \*Kosa, G., Shoham, M. and Zaaroor M., IEEE Conference on Robotics and Automation (ICRA05), (2005) 1339-1343.

Brigham and Women's Hospital, Boston


Nobuhiko Hata

## Our approach

- To combine a magnetic coil and swimming tails and convert the force-couple created by the magnetic coil into a propulsive force in the direction of the static magnetic field
- We will design and develop an optimal waving pattern for the tails and try to achieve advancing speed of 5 mm/sec with a magnetic force of 5 mN.
- The capsule will be 3 mm in diameter and 20 mm long including the three swimming tails.

Brigham and Women's Hospital, Boston


Nobuhiko Hata



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Prototype





Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Simulation

- Length of swimming tail 10 [mm]
- Width 5 [mm]
- Total thickness 0.15 [mm]
- Current 3.7 [mA] in the MRI's magnetic field 3 [T].

- Traveling wave
- Amplitude 14 [mm]
- Frequency of 12.836 [kHz]



Brigham and Women's Hospital, Boston


Nobuhiko Hata

## Future direction

- Miniaturization by MEMS
- Treatment
- Molecular imaging

Brigham and Women's Hospital, Boston

Nobuhiko Hata




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Engineering

Brigham and Women's Hospital, Boston

Nobuhiko Hata




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## What is 3D Slicer?

- 3D Slicer is...
  - An end-user application for 3D medical image computing research and Image Guided Therapy
  - A platform for research where new techniques can be plugged into a useful framework
  - A freely-downloadable program with source and binaries for Windows, Linux, Solaris and (increasingly) Mac OS



Slide courtesy of Drs. Kikinis and Pieper

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Slicer Background

- SPL Image Guided Surgery and Visualization (Kikinis, Westin, Hata, Halle, others)
- Slicer Application Pulled Together by Dave Gering 1997-1999 with VTK and Tcl
- Further Development and Architecture by Lauren O'Donnell 1999-2001
- Ongoing Development of Slicer's Base Primarily by Steve Pieper and Nicole Aucoin
- **Many Modules and Contributions by Various Authors**
  - BWH, MIT, MGH, Georgia Tech, UCSD, JHU...
- Now being used as a reference application platform for NA-MIC

•Slide courtesy of Drs. Kikinis and Pieper

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Features

- **Load Medical Image Data: MR, CT in DICOM, GE, Analyze...**
- **XML-based File Format: MRML (Medical Reality Markup Language)**
- **Interactive Editor: Draw, Threshold, Math Morphology...**
- **Automated Segmenters: EM Segmentation, Fast Marching, Level Sets...**
- **Visualization: Model Building, Stereo Rendering, Animation...**
- **Registration: Manual, ITK, CNI**
- **Measurement: Fiducial-Based, Volumetric, Polyhedral Intersection, Vessel Cross-Section, Osteotomy Planning**
- **IGT: Tracked Probes, Real-Time Images, Robot Control**
- **OPEN SOURCE**


•Slide courtesy of Drs. Kikinis and Pieper

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Slicer in IGT

- To facilitate application of state-of-art medical image processing in IGT
- IGT specific extension to support this translation
  - Tracker (IGSTK, Cleary RO1)
  - GE MRI scanner
  - Endoscope support
  - (3D US)
  - (Siemens, Philips)




Brigham and Women's Hospital, Boston

Nobuhiko Hata

## IGT Software Design

- Software design to maximize function commonalities among applications
  - Brain (biopsy, craniotomy, NdYAG laser ablation)
  - Prostate (brachytherapy, biopsy)
  - Liver and kidney (Microwave, Cryo, laser ablation)
  - Endoscopy (broncho-, neuro-, feto-scopy)



Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Open-source software: why?

- **Resource Sharing**
  - Lower development cost
  - Dissemination and training
  - Facilitates communication and collaboration
- **Open architecture**
  - Access to proprietary hardware
- BSD-style license for potential commercialization

Brigham and Women's Hospital, Boston

Nobuhiko Hata

## Conclusion

- Robot, as positioning tool, enhances the precision of the MRI-guided therapy
  - Autonomous -> Interactive
- Robot can be an enabling tool to perform therapies in close-bore scanner
- Added value - organ motion compensation will enable IGT of moving organ
- Capsule endoscope, with an addition of image guidance, can be extended to treatment

Brigham and Women's Hospital, Boston

Nobuhiko Hata

- Slides:
  - [www.snrlab.org](http://www.snrlab.org) -> events